



CLVision™

Angiography Catheter

CLVision™

Clear Visualization. Exceptional Control.

CLVision™ Angiography Catheter is engineered to deliver high-quality contrast media injection while providing excellent torque control and vessel access. The reinforced braided shaft combines flexibility with pushability, making it suitable for a wide range of diagnostic angiography procedures.

Key Features

- Pebax shaft with stainless-steel braided reinforcement.
- Outstanding pushability and trackability.
- Designed to achieve 1:1 torque.
- Soft atraumatic tip minimizes vessel trauma.
- Large inner lumen for high contrast flow rates.
- Pressure resistance up to 1000 psi.
- Radiopaque marker for enhanced fluoroscopic visibility for some variants.
- Multiple curve options to support physician preference.

TECHNICAL OVERVIEW

French Sizes	4F • 5F • 6F • 7F
Working Length	100 cm, 110 cm
Tip Design	Soft Atraumatic Tip
Maximum guidewire size:	0.038"/0.97mm
Visibility	Radiopaque Markers (Pig tail)
Usage	Single Use

Coronary Engagement Catheters

Judkins Left (JL)
Judkins Right (JR)
Amplatz Left (AL)
Amplatz Right (AR)

Specialty Catheters


Coronary Bypass Left (CBL)
Coronary Bypass Right (CBR)
Internal Mammary (IM)
Tiger (TIG)

Multipurpose Catheters

MPA
MPB

Ventriculography & Aortography Catheters

Pigtail (PIG)



Ordering Information

Judkins Left	Catheter Length	Catheter Size	Catalogue Number
JL3.0	100cm	4F	CLJL30I004F
JL3.0	100cm	5F	CLJL30I005F
JL3.0	100cm	6F	CLJL30I006F
JL3.0	100cm	7F	CLJL30I007F
JL3.5	100cm	4F	CLJL35I004F
JL3.5	100cm	5F	CLJL35I005F
JL3.5	100cm	6F	CLJL35I006F
JL3.5	100cm	7F	CLJL35I007F
JL4.0	100cm	4F	CLJL40I004F
JL4.0	100cm	5F	CLJL40I005F
JL4.0	100cm	6F	CLJL40I006F
JL4.0	100cm	7F	CLJL40I007F
JL4.5	100cm	4F	CLJL45I004F
JL4.5	100cm	5F	CLJL45I005F
JL4.5	100cm	6F	CLJL45I006F
JL4.5	100cm	7F	CLJL45I007F
JL5.0	100cm	4F	CLJL50I004F
JL5.0	100cm	5F	CLJL50I005F
JL5.0	100cm	6F	CLJL50I006F
JL5.0	100cm	7F	CLJL50I007F
JL6.0	100cm	4F	CLJL60I004F
JL6.0	100cm	5F	CLJL60I005F
JL6.0	100cm	6F	CLJL60I006F
JL6.0	100cm	7F	CLJL60I007F

Judkins Right	Catheter Length	Catheter Size	Catalogue Number
JR3.0	100cm	4F	CLJR30I004F
JR3.0	100cm	5F	CLJR30I005F
JR3.0	100cm	6F	CLJR30I006F
JR3.0	100cm	7F	CLJR30I007F
JR3.5	100cm	4F	CLJR35I004F
JR3.5	100cm	5F	CLJR35I005F
JR3.5	100cm	6F	CLJR35I006F
JR3.5	100cm	7F	CLJR35I007F
JR4.0	100cm	4F	CLJR40I004F
JR4.0	100cm	5F	CLJR40I005F
JR4.0	100cm	6F	CLJR40I006F
JR4.0	100cm	7F	CLJR40I007F
JR4.5	100cm	4F	CLJR45I004F
JR4.5	100cm	5F	CLJR45I005F
JR4.5	100cm	6F	CLJR45I006F
JR4.5	100cm	7F	CLJR45I007F
JR5.0	100cm	4F	CLJR50I004F
JR5.0	100cm	5F	CLJR50I005F
JR5.0	100cm	6F	CLJR50I006F
JR5.0	100cm	7F	CLJR50I007F
JR6.0	100cm	4F	CLJR60I004F
JR6.0	100cm	5F	CLJR60I005F
JR6.0	100cm	6F	CLJR60I006F
JR6.0	100cm	7F	CLJR60I007F

Amplatz Left	Catheter Length	Catheter Size	Catalogue Number
AL1	100cm	4F	CLAL1004F
AL1	100cm	5F	CLAL1005F
AL1	100cm	6F	CLAL1006F
AL1	100cm	7F	CLAL1007F
AL2	100cm	4F	CLAL21004F
AL2	100cm	5F	CLAL21005F
AL2	100cm	6F	CLAL21006F
AL2	100cm	7F	CLAL21007F
AL3	100cm	4F	CLAL31004F
AL3	100cm	5F	CLAL31005F
AL3	100cm	6F	CLAL31006F
AL3	100cm	7F	CLAL31007F
AL4	100cm	4F	CLAL41004F
AL4	100cm	5F	CLAL41005F
AL4	100cm	6F	CLAL41006F
AL4	100cm	7F	CLAL41007F

Amplatz Right	Catheter Length	Catheter Size	Catalogue Number
AR1	100cm	4F	CLAR11004F
AR1	100cm	5F	CLAR11005F
AR1	100cm	6F	CLAR11006F
AR1	100cm	7F	CLAR11007F
AR2	100cm	4F	CLAR21004F
AR2	100cm	5F	CLAR21005F
AR2	100cm	6F	CLAR21006F
AR2	100cm	7F	CLAR21007F
AR3	100cm	4F	CLAR31004F
AR3	100cm	5F	CLAR31005F
AR3	100cm	6F	CLAR31006F
AR3	100cm	7F	CLAR31007F

CBR	Catheter Length	Catheter Size	Catalogue Number
CBR	100cm	4F	CLCBRI004F
CBR	100cm	5F	CLCBRI005F
CBR	100cm	6F	CLCBRI006F
CBR	100cm	7F	CLCBRI007F

CBL	Catheter Length	Catheter Size	Catalogue Number
CBL	100cm	4F	CLCBLI004F
CBL	100cm	5F	CLCBLI005F
CBL	100cm	6F	CLCBLI006F
CBL	100cm	7F	CLCBLI007F

IM	Catheter Length	Catheter Size	Catalogue Number
IM	100cm	4F	CLIMI004F
IM	100cm	5F	CLIMI005F
IM	100cm	6F	CLIMI006F
IM	100cm	7F	CLIMI007F

MPA	Catheter Length	Catheter Size	Side Holes	Catalogue Number
MPA1	100cm	4F	0	CLMPA11004F
MPA1	100cm	5F	0	CLMPA11005F
MPA1	100cm	6F	0	CLMPA11006F
MPA1	100cm	7F	0	CLMPA11007F
MPA2	100cm	4F	2	CLMPA21004F
MPA2	100cm	5F	2	CLMPA21005F
MPA2	100cm	6F	2	CLMPA21006F
MPA2	100cm	7F	2	CLMPA21007F

MPB	Catheter Length	Catheter Size	Side Holes	Catalogue Number
MPB1	100cm	4F	0	CLMPB11004F
MPB1	100cm	5F	0	CLMPB11005F
MPB1	100cm	6F	0	CLMPB11006F
MPB1	100cm	7F	0	CLMPB11007F
MPB2	100cm	4F	2	CLMPB21004F
MPB2	100cm	5F	2	CLMPB21005F
MPB2	100cm	6F	2	CLMPB21006F
MPB2	100cm	7F	2	CLMPB21007F

PIG	Catheter Length	Catheter Size	Side Holes	Catalogue Number
PIG	110cm	4F	6	CLPIG1104F
PIG	110cm	5F	6	CLPIG1105F
PIG	110cm	6F	5	CLPIG1106F
PIG	110cm	7F	5	CLPIG1107F
PIG A 145°	110cm	5F	6	CLPIGA1105F
PIG A 145°	110cm	6F	5	CLPIGA1106F
PIG B 155°	110cm	5F	6	CLPIGB1105F
PIG B 155°	110cm	6F	5	CLPIGB1106F

TIG	Catheter Length	Catheter Size	Side Holes	Catalogue Number
TIG	100cm	4F	1	CLTIG1004F
TIG	100cm	5F	1	CLTIG1005F
TIG	100cm	6F	1	CLTIG1006F
TIG	100cm	7F	1	CLTIG1007F



MedivationBio

Achieving Possibilities

Medivation Bio is an emerging medical device company committed to delivering safe, reliable and high-quality consumables for critical care and routine medical use. Our product portfolio includes IV Cannulas, Hypodermic Syringes, Blood Collection Tubes, Chest Drainage Catheters, High-Pressure Tubings and Cardiovascular Devices.

We follow strict quality protocols and design products on:

Patient safety

Reliable engineering standards

Ease of use

Clinical performance

"Achieving Possibilities"

is our commitment to creating devices that support clinicians and improve patient care.

SIDCO Industrial Complex, Bari Brahmana EPIP Industrial Estate,
Samba, Jammu and Kashmir, India 181133.