



# RazorTip™

# Introducer Sheath


## Precision Access Starts Here

Designed to facilitate safe and efficient vascular access, the RazorTip™ Introducer Set provides smooth catheter introduction, reliable hemostasis, and excellent procedural control across a broad range of diagnostic and interventional applications.

- Atraumatic tapered Sheath and dilator tip for smooth vascular access.
- Advanced cross-cut hemostatic valve technology ensures effective hemostasis, supporting procedural efficiency and operator confidence.
- Refined sheath-dilator interface for reduced insertion resistance and enhanced ease of access.
- Colour-coded hub enables instant identification of sheath French size during procedures.
- Side port with 3-way stopcock enables convenient flushing, pressure monitoring and contrast administration.
- Secure locking mechanism ensures stable sheath-dilator engagement during insertion.
- Available in femoral and transradial configurations.
- Wide range of French sizes, sheath lengths and guidewire combinations to meet diverse clinical requirements.

## Technical Specifications

Parameter	Available Options
Sheath Size	5F, 6F, 7F, 8F
Sheath Length	7 cm, 11 cm
Guidewire Options	0.018", 0.021", 0.035", 0.038"
Guidewire Type	Straight Tip / J-Tip
Needle Options	18G, 21G
Sterility	EO Sterilized



## RazorTip™ Introducer set

Introducer Access Type	Sheath Size	Sheath Length	Dilator	Guidewire O.D. x Length	Needle size	Catalogue Number	
						S tip	J tip
Femoral	5F	11cm	16cm	0.035" x 45cm	18G x 7cm	RT5F11351807S	RT5F11351807J
	6F					RT6F11351807S	RT6F11351807J
	7F					RT7F11351807S	RT7F11351807J
	8F					RT8F11351807S	RT8F11351807J

Introducer Access Type	Sheath Size	Sheath Length	Dilator	Guidewire O.D. x Length	Needle size	Catalogue Number	
						S tip	J tip
Femoral	5F	11cm	16cm	0.038" x 45cm	18G x 7cm	RT5F11381807S	RT5F11381807J
	6F					RT6F11381807S	RT6F11381807J
	7F					RT7F11381807S	RT7F11381807J
	8F					RT8F11381807S	RT8F11381807J

Introducer Access Type	Sheath Size	Sheath Length	Dilator	Guidewire O.D. x Length	Needle size	Catalogue Number	
						S tip	J tip
Transradial	5F	11cm	16cm	0.021" x 45cm	21G x 4cm	RT5R11212104S	RT5R11212104J
	6F			RT6R11212104S		RT6R11212104J	
	5F			0.018" x 45cm		RT5R11182104S	RT5R11182104J
	6F			RT6R11182104S		RT6R11182104J	

Introducer Access Type	Sheath Size	Sheath Length	Dilator	Guidewire O.D. x Length	Needle size	Catalogue Number	
						S tip	J tip
Transradial	5F	7cm	13cm	0.021" x 45cm	21G x 4cm	RT5R07212104S	RT5R07212104J
	6F			RT6R07212104S		RT6R07212104J	
	5F			0.018" x 45cm		RT5R07182104S	RT5R07182104J
	6F			RT6R07182104S		RT6R07182104J	



# MedivationBio

Achieving Possibilities

Medivation Bio is an emerging medical device company committed to delivering safe, reliable and high-quality consumables for critical care and routine medical use. Our product portfolio includes IV Cannulas, Hypodermic Syringes, Blood Collection Tubes, Chest Drainage Catheters, High-Pressure Tubings and Cardiovascular Devices.

We follow strict quality protocols and design products on:

Patient safety

Reliable engineering standards

Ease of use

Clinical performance

"Achieving Possibilities"

is our commitment to creating devices that support clinicians and improve patient care.

SIDCO Industrial Complex, Bari Brahmana EPIP Industrial Estate,  
Samba, Jammu and Kashmir, India 181133.